



Appearance:

Stream-line design

Water drainage design

Compatible with other standard chambers

流線、排水設計，並相容於其他加濕瓶

HEMOCARE HUMIDIFIER / 居家型加熱器

REF / 型號	VH-2000-115 / VH-2000-230
Dimensions / 規格	W136 x D160 x H61 mm
Weight / 重量	0.9 kg (w/o Chamber fitted)
Electrical Power / 電源	(1) Heater plate capacity : 150W max (2) Temperature Power cutout at 115°C ± 3°C (3) 230 ± 23V, 0.7A max. at 230V (4) 115 ± 12V, 1.3 A max. at 115V
Heater Plate Overheat / 過熱保護	Cutoff: 93±5°C
Temperature Setting / 溫度設定	(1) From step 1 to step 5 (2) Range 45~65°C
Standards / 符合標準	(1) IEC/EN 60601-1, IEC/EN 60601-1-2 (2) ISO 8185
Accessories / 配件	VH-3147 Reu. Humidifier Chamber w/o Scroll
Package / 包裝尺寸	12 sets per carton Carton: L41 x W41 x H29cm