

CLEARVUE™

VIDEO LARYNGOSCOPE



INFINIUM 

 AI-KARAM
Medical Devices

CLEARVUE™



The Infinium ClearVue™ VL3R

Reusable Blade Video Laryngoscopes

The Infinium ClearVue™ video laryngoscope adopts state-of-the-art image technology and possesses the advantage of a 2.0 megapixel full view camera with 4 inch high-resolution monitor. The ClearVue™ features a rechargeable li-ion battery for extended use. The HD screen offers a unique anti-fog and anti-shatter capability. The overall design is compact and ergonomic. A high-resolution view of the glottis will enable clinicians first-attempt success while minimizing any chance of a complication during the intubation process.

The Infinium ClearVue™
Video Laryngoscopes are:

COMPACT AND ERGONOMIC

- Lightweight - Unit Is Less Than 350g
- Rigid Steel Construction
- Quick installation and Removal of Blades
- Quick Shot Camera Button for Photo and Video Recording
- 4 GB Data Storage (Approx. 80 Minutes of Stored Video)
- Optional HDMI Output for Remote Display

COMPREHENSIVE BLADE OFFERING

- Macintosh Blades Size 1,2,3,4,5 Available
- Macintosh Curved “D” Blades Size 1,2,3,4,5 Available
- Miller Blades Size 0,1 available
- The ClearVue Comes standard with Macintosh Size 2,3, and 4 Blades




DURABLE

- 2 Years Warranty
- Screen Reinforced with Anti-Shatter Protection
- Light Weight Long Lasting Rechargeable Li-ion Battery

HIGH-DEFINITION

- 4 Inch HD Display
- Anti-Fog Lens
- 2.0 Mega Pixel Camera With Video

Reusable Blade Attachment Options



The ClearVue™ VL3R comes standard with size 2, 3, and 4 reusable Mac blades.



Optional blade sizes available are:
Miller 0



Miller 1



Mac 1



and Mac 5.



Optional Curved Mac "D" Blades
Sizes 1,2,3,4,5

CLEARVUE® TECHNICAL SPECIFICATIONS:

DISPLAY

Size	4" Full View
Monitor	LCD; Scratch-resistant and anti-glare screen
Resolution Ratio	640*480 RGB
Aspect Ratio	4:3
Video Refresh Rate	30FPS
Illumination	LED

LARYNGOSCOPE BLADE CAMERA INCLUDED

Resolution Ratio	2.0M Pixels
Field Angle	66°
Illuminance	≥800LUX

BATTERY

Type	Rechargeable Lithium Battery
------	------------------------------

Long working time	≥120minutes
Voltage	3.7V
Capacity	3200mAh
Charging times	More than 300 times
Charging time	4 hours

POWER ADAPTER

Charging port	Micro USB
Input	100-250V,50/60Hz
Output	5V/ 2A

TRANSPORT/ STORAGE

Temperature	-10 ~+45
-------------	----------

CONDITION

Humidity	≤93
AP	500hPa to 1060hPa

WORKING CONDITION

Temperature	5 ~+40
Humidity	30 ~85
AP	700hPa to 1060hpa

UNIQUE ANTI-FOG FUNCTION

Auto anti-fog function upon powering on

WEIGHT

321 g(disposable system); 376g (reusable system)
- light-weight, portable, comfortable

PHOTO/VIDEO FUNCTION

One click image capture 4G internal memory card for video and images USB connection cable for uploading images and video

HANDLE

ABS material - Comfortable ergonomic short handle design

CLEARVUE ORDER INFORMATION & PART NUMBERS

000.300.0000	Laryngoscope kit with a set of 3 reusable blades (VL3R)	261.300.0004	Laryngoscope reusable blade No.4
000.300.0002	Laryngoscope kit with a set of 3 reusable blades and HDMI cable (VL3RH)	261.300.0005	Laryngoscope reusable blade No.5
000.300.0003	Laryngoscope kit with a set of disposable blades (new version)	261.300.0006	HDMI to microUSB cable kit(with adapter) for laryngoscope
000.300.0004	ClearVue III, 8-in display Laryngoscope kit with a set of disposable blades	261.300.0007	ClearVue Reusable Blade Difficult Airway
000.300.0005	ClearVue Reusable Handle (Handle Only, No Blades)	261.300.0008	ClearVue Reusable Blade MILLER Size 0
000.300.0006	Laryngoscope 4in display kit with a set of 3 reusable blades	261.300.0009	ClearVue Reusable Blade MILLER Size 00
000.300.0007	Laryngoscope 4in display kit with a set of 3 reusable blades and HDMI cable	261.300.0010	ClearVue Reusable Blade MILLER Size 1
000.300.0008	Laryngoscope 4in display kit with a set of disposable blades (new version)	261.300.0011	Reusable D Blade (60) Mac 1
		261.300.0012	Reusable D Blade (60) Mac 2
261.300.0001	Laryngoscope reusable blade No.1	261.300.0013	Reusable D Blade (60) Mac 3
261.300.0002	Laryngoscope reusable blade No.2	261.300.0014	Reusable D Blade (60) Mac 4
261.300.0003	Laryngoscope reusable blade No.3	261.300.0015	Reusable D Blade (60) Mac 5



10th of Ramadan City, Alordonia, Al Jawhara Mall, Office 447.
32 Obur Bulding - Salah Salem Street, Cairo, Egypt.

sales@alkaram-medical.com
www.alkaram-medical.com
01099442054 - 01206467528
01272398884 - 01032671638



12151 62nd St N #5 • Largo FL 33773 • USA
Phone: (1) 727-531-8434 • Fax: (1) 727-531-8436
Web: www.infiniummedical.com • Email: sales@infiniummedical.com

043.300.0002 Rev 1